



Titan Chlor Plus Tabs


2 in 1 Cleaner and Disinfectant, Chlorine Releasing Tablets (1.7g Troclosene)

Description

DTitan Chlor Plus Tabs are a blend of NaDCC, effervescent agents and tableting agents. They provide an accurate and convenient means of making up solutions of known available chlorine strengths for a variety of disinfection tasks. Titan Chlor PlusTabs are soluble in warm or cold water and are compatible when dissolved with neutral detergents for making a one stage sanitising solution. Contains 1.7g Troclosene

Use instructions

Environmental disinfection: For increased safety in general cleaning of high risk areas (theatres, maternity units, etc.) drains and sinks, use 1 tablet per litre of water of neutral detergent (1000 ppm available chlorine). Refer to dilution table for more information.



© A.I.S.E.



Titan Chlor Plus Tabs

Technical data

Appearance: firm white tablet

Odour: Neutral

pH value (in use): 1.95

Relative density (20°C): 1.000

Flash point (°C): N/A

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diverseys.com. Store in original closed containers away from extremes of temperature. Only for professional users / specialists.

12 months minimum under recommended storage conditions

Microbiological data

Chlorine release agents have been shown to be effective against a wide range of micro-organisms: bacteria, yeast, fungi, viruses, spores.

Efficacy confirmed against European norms: EN 1276, EN 13727, EN 13697, EN 14476
EN1650, EN13624, EN13704, EN17126

Available pack sizes

Titan Chlor PlusTabs are available in tubs of 200 tablets, 6 tubs per case.

Medguard Professional Healthcare

Tel: 01 835 2378

Email: orders@medguard.ie

Fax: 01 969 5009

Web: www.medguard.ie

v

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, Ireland



Titan Chlor Plus Tabs

1 tablet in	Gives Available Chlorine of	Typical uses
100 mls	10,000 ppm	Body fluid spils and conditions of heavy soiling
400 mls	2,500 ppm	Non metallic medical and Vetinary appliances and Laboratory discard jars
1 litre	1,000 ppm	General disinfection
2.5 litres	400 ppm	W.C's, drains, sinks
5 litres	200 ppm	Food preperation Areas
10 litres	100 ppm	Cloths and mops