

# seca 203 & seca 201

## Ergonomic Measuring Tape

new



- Retractable measuring tape for precise body measurements.
- WHR calculator for quick determination of body fat distribution (only 203).
- WHR scale of metal for long life (only 203).

Medguard Healthcare

Tel: +353 (01)835 24 11 Email: [service@medguard.ie](mailto:service@medguard.ie)

Fax: +353 (01)835 24 61 Web: [www.medguard.ie](http://www.medguard.ie)

## seca 203 & seca 201: Two high-quality medical products.

In the seca 203 and 201 you've found exactly the right products if you're in search of a measuring tape that is clearly a cut above the rest. Or if you're looking for a measuring tape with millimeter markings that can accurately measure to five millimeters.

### seca 203:

The seca 203 also has an extra instrument to calculate the Waist-to-Hip Ratio (WHR). This ratio is an indirect method of determining the amount of fatty tissue in the abdomen and consequently, the distribution of fat reserves. The high-quality scale is made of metal so that you'll be able to depend on this valuable additional function for a long time.



### seca 201:

The seca 201 is a true medical product. It's very easy to pull out the tape, which just as easily clicks back into place. Handling is extraordinarily convenient. Made of durable



materials, the scale features a tape with a top measurement of 205 centimeters and its own low weight of only 55 grams.



seca 203:  
Display to determine  
body fat distribution.



seca 201:  
High-quality  
materials for precise  
functionality.

## seca 203 & seca 201

### Technical Data

#### seca 203:

- Measurement range: 15 - 205 cm
- Division: 1mm
- Diameter: 70 mm
- Height: 28 mm
- Weight of device: about 75 g
- Function: WHR calculator

#### seca 201:

- Measurement range: 15 - 205 cm
- Division: 1 mm
- Diameter: 70 mm
- Height: 22 mm
- Weight of device: about 55 g

Medguard Healthcare

Tel: +353 (01)835 24 11 Email: [service@medguard.ie](mailto:service@medguard.ie)

Fax: +353 (01)835 24 61 Web: [www.medguard.ie](http://www.medguard.ie)