

seca 799

Digital column scale
with BMI function 



- Site-independent as battery operated.
- Easy to transport.
- Cost-effective and low-maintenance due to low power consumption.
- Intelligent functions for versatile use.

Medguard Healthcare

Tel: +353 (01)835 24 11 Email: service@medguard.ie

Fax: +353 (01)835 24 61 Web: www.medguard.ie

seca 799:

The intelligent solution for flexible use

A must in the medical sector: high precision and sophisticated functionality combined with cost-effectiveness. For more than one hundred years now, seca has been developing innovative scales and measuring systems to provide optimal solutions. Because only equipment which is perfectly aligned to the requirements of doctors and nursing staff and which, at the same time, work cost-effectively can really ease the workload in hospitals, practices and nursing homes. The versatile column scale seca 799 is one of these developments.



Intelligent functions for versatile use.

Essential in hospitals and doctors' practices: the determination of the nutritional condition of the patient. That is why the seca 799 has the BMI function where the height of the patient is compared with his weight, thus providing valuable information. The HOLD function too reduces the everyday workload: it ensures that the result is still shown on the display after the patient has stepped down from the scale. The carer can therefore first attend to the patient before making a note of the weight. Additional items which are not to be included in the weight result are simply deducted by means of the TARE function

Site-independent and easy to transport.

The seca 799 can be used wherever it is needed. As it is battery-operated, no wall outlet connection is necessary. And a versatile scale must also be easy to move around – that is why the seca 799 is fitted with smooth running transport castors on which it can quite simply be wheeled to wherever it is needed.

Especially cost-effective and low-maintenance.

The seca 799 is extremely energy-saving and thus very cost-effective. The minimum power consumption of the low maintenance column scale ensures that a large number of weighings are possible with just one set of batteries. The automatic switch-off function is another reason for the low energy consumption.



seca 799

Medguard Healthcare

Tel: +353 (01)835 24 11 Email: service@medguard.ie

Fax: +353 (01)835 24 61 Web: www.medguard.ie

Technical Data

- Capacity: 200kg
- Graduation: 100 g < 150 kg > 200 g
- Dimensions (WxHxD):
294 x 831 x 417 mm
- Dimensions/platform (WxHxD):
272 x 75 x 280 mm
- Weight: 6.3 kg
- Power supply: batteries, mains adapter optional
- Functions: TARE, HOLD, BMI, weighing range switch-over, automatic switch-off
- Optional: measuring rod seca 220, seca 224, mains adapter seca 400

Medguard Healthcare

Tel: +353 (01)835 24 11 Email: service@medguard.ie

Fax: +353 (01)835 24 61 Web: www.medguard.ie