



Stand scales

for pharmacies and care homes

7830

Stand scale

SOEHNLE
PROFESSIONAL

- **Comfortable display height**
- **Improved platform surface**

- Bright luminous LED display, 15 mm digit height
- Extra large tread
- Instant-start-function by stepping on the scale
- Tare function: Button on the back of the glass display
- Mother-child-weighing
- Very easy to clean surface
- Interface for PC/IT usable with data cable 2550.03.019
- Units in kg or lb
- Energy-saving auto-off function
- Mains operation



Order number	7830.30.001
Maximum capacity	250 kg
Division	100 g
Dimensions	365 x 490 x 890 mm
Scale weight	14,1 kg

Accessories

Data cable serial and USB interface for data transfer PC/IT
Order number 2550.03.019

Rechargeable battery, premounted*
Order number 2557.12.002

Rolls (2 pieces), premounted*
Order number 2513.12.002

*The installation of the rechargeable battery and rolls can be carried out only at the factory.



7831

Stand scale with height rod

SOEHNLE
PROFESSIONAL

S *Automatic BMI calculation*

S *Improved platform surface*

- > Integrated height rod
- > Indication of weight, height and automatic BMI calculation
- > Bright luminous LED display, 15 mm digit height
- > Instant-start-function by stepping on the scale
- > Tare function: Button on the back of the glass display
- > Mother-child-weighing
- > Very easy to clean surface
- > Interface for PC/IT usable with data cable 2550.03.019
- > Units in kg or lb
- > Energy-saving auto-off function
- > Mains operation



Order number	7831.30.001
Maximum capacity	250 kg
Division	100 g
Measuring range	120 - 202 cm
Dimensions	365x490x1,280 mm
Scale weight	17.2 kg

Accessories

Data cable serial and USB interface for data transfer PC/IT

Order number 2550.03.019

Rechargeable battery, premounted*

Order number 2557.12.002

Rolls (2 pieces), premounted*

Order number 2513.12.002

* The installation of the rechargeable battery and rolls can be carried out only at the factory.

Medguard

Professional Healthcare Supplies

Tel: +353 (01) 835 2411 Email: service@medguard.ie

Fax: +353 (01) 835 2461 Web: www.medguard.ie

