



7730

SOEHNLE
PROFESSIONAL

Stand scale approved

S *Comfortable display height*

S *Improved platform surface*

- › Approved to class III, Medical device directive (MPG) class 1
- › Extra large tread
- › Bright luminous LED display, 15 mm digit height
- › Instant-start-function by stepping on the scale
- › Tare function: Button on the back of the glass display
- › Mother-child-weighing
- › Very easy to clean surface
- › Interface for PC/IT usable with data cable 2550.03.019 as service interface
- › Units in kg or lb
- › Energy-saving auto-off function
- › Mains operation



Order number	7730.01.001
Maximum capacity	250 kg
Division	50 g / 100 g
Dimensions	365x490x890 mm
Scale weight	14,1 kg

Accessories

Data cable for service interface
Order number 2550.03.019

Rechargeable battery, premounted*
Order number 2557.12.002

Rolls (2 pieces), premounted*
Order number 2513.12.002

* The installation of the rechargeable battery and rolls can be carried out only at the factory.

Approved Stand scales for hospitals and doctors



7731

Stand scale with height rod approved

SOEHNLE
PROFESSIONAL

S Automatic BMI calculation

S Improved platform surface

- › Approved to class III, Medical device directive (MPG) class 1
- › Integrated height rod
- › Indication of weight, height and automatic BMI calculation
- › Bright luminous LED display, 15 mm digit height
- › Instant-start-function by stepping on the scale
- › Tare function: Button on the back of the glass display
- › Mother-child-weighing
- › Very easy to clean surface
- › Interface for PC/IT usable with data cable 2550.03.019 as service interface
- › Units in kg or lb
- › Energy-saving auto-off function
- › Mains operation



Order number	7731.01.001
Maximum capacity	250 kg
Division	50 g / 100 g
Measuring range	120 - 202 cm
Dimensions	365x490x1,280 mm
Scale weight	17.2 kg

Accessories

Data cable for service interface
Order number 2550.03.019

Rechargeable battery, premounted*
Order number 2557.12.002

Rolls (2 pieces), premounted*
Order number 2513.12.002

* The installation of the rechargeable battery and rolls can be carried out only at the factory.

Medguard
Professional Healthcare Supplies
Tel: +353 (01) 835 2411 Email: service@medguard.ie
Fax: +353 (01) 835 2461 Web: www.medguard.ie

