



7732

SOEHNLE
PROFESSIONAL

MediScale approved

- S** Various interfaces for transfer of weight values to a PC, printer or to a mobile phone
- S** Improved platform surface

- > Extra large tread
- > Approved to class III, Medical device directive (MPG) class 1
- > Bright luminous LED display, 15 mm digit height
- > Instant-start-function by stepping on the scale
- > Units in kg or lb
- > Energy-saving auto-off function
- > Mains operation



Personal scales for rehabilitation clinics and sport centres

Order number	7732.01.001	7732.30.001
Model	Standard	RS232
		Data transfer
Maximum capacity	250 kg	250 kg
Division	50 g / 100 g	50 g / 100 g
Dimensions	365x490x170 mm	365x490x170 mm
Scale weight	13 kg	13 kg
Order number	7732.10.001	OEM
Model	Bluetooth™	OEM
	Wireless transmission	Custom Model
Maximum capacity	250 kg	250 kg
Division	50 g / 100 g	50 g / 100 g
Dimensions	365x490x170 mm	365x490x170 mm
Scale weight	13 kg	13 kg

Medguard

Professional Healthcare Supplies

Tel: +353 (01) 835 2411 Email: service@medguard.ie

Fax: +353 (01) 835 2461 Web: www.medguard.ie

