



Chair weighers

for hospitals and nursing homes




Basic version: Display Basic / Basic E




Optional: Display Standard approveable



Optional: Display Glass



By folding up the armrests, getting onto the chair from the hospital bed is made easier.



Hinged footrests for easy entry and exit.

7702/7802

Chair weigher

S Robust construction

S Simple weighing of persons with a mobility handicap

- › Tare function: Covers/paddings can be tared
- › Hold function for better readability
- › Ergonomically shaped and easy-to-clean seat
- › Hinged arm- and footrests
- › 2 wheels on backside are rotatable and with parking brake
- › Units in kg or lb
- › Energy-saving auto-off function
- › Rechargeable battery operation with integrated charger (operating time 50 h, charging time 5.5 h)
- › Display freely selectable: Basic, Glass, Basic E or Standard (for a detailed description of the displays, see datasheet "Displays")



7702

Chair weigher approved

Approval	approved to class III, MPG class 1
Maximum capacity	200 kg
Division	50 g / 100 g
Dimensions	580x1,100x930 mm
Scale weight	25 kg

With display Basic E

Order number 7702.37.002

With display Standard

Order number 7702.65.002



7802

Chair weigher

Maximum capacity	200 kg
Division	100 g
Dimensions	580x1,100x930 mm
Scale weight	25 kg

With display Basic

Order number 7802.70.002

With display Glass

Order number 7872.01.002



Accessories

Bag for charger

Order number 2507.03.001

Medguard
Professional Healthcare Supplies
Tel: +353 (01) 835 2411 Email: service@medguard.ie
Fax: +353 (01) 835 2461 Web: www.medguard.ie

Interface and data cable for transfer of weight values, see datasheet "Displays"