


- Temp range: +2°C/+8°C
- External digital controller accurate to 0.1°C
- Forced air cooling
- Electronic controller
- Max & min temperature recording
- High & low alarm
- Visual & Audible alarm
- Glass door
- Lockable

### Special Features

- Safety thermostat to prevent freezing of medicaments
- Door alarm
- LED low energy internal lighting
- Auxiliary output for remote alarm / BMS system



### Option

- Solid Door

ARVF 300

### Options

- Drawer with dividers

Model	Dimensions ( L x D x H mm )	Capacity ( L )	Shelves
ARVF 300	595 x 608 x 1563	278	6
Cool Tray Drawer*	457 x 404 x 60	-	-

\*with 3 lanes & 12 internal divisions - Out of 2 years warranty