



## Dopplers

Dopplers are available as Fetal or Vascular. Fetal detect the fetal heart rate (FHR) and Vascular determine the amount of blood flow through a vessel.

Most dopplers are for obstetric use (fetal monitoring) and come with a fixed probe. This means that you cannot change the probe at anytime.

Some dopplers allow interchangeable probes so you can alternate between a fetal or a vascular probe.

Note that some dopplers come with a probe and some have to be purchased seperately.

Fetal Dopplers usually come with a 2Mhz or 3Mhz probe:-

- 2Mhz - More for general use and commonly detect FHR around 10 weeks.
- 3Mhz - More enhanced, offering detection of FHR around 8 weeks.

Vascular Doppler Probes are available in the following Mhz:-

- 4Mhz - Detection of deep vessels.
- 5Mhz - Commonly used for vascular testing of deep lying vessels above the knee or elbow.
- 8Mhz - Most popular vascular probe, used for peripheral vessels below the knee and below the elbow.
- 10Mhz - Specialist probe commonly used by plastic surgeons.



### Sonotrax Doppler

Can be used as Fetal Doppler or Vascular Doppler. Comes with choice of four probes, 2,3,4 or 8 Mhz. Basic audio quality. Visual display of FHR.



### Pocket Fetal Doppler

Fetal doppler with 2Mhz probe. Probe can be changed but we are promoting with 2Mhz only. Only suitable after 12th week of pregnancy. Audio and visual display of FHR. Basic quality.



### Sonicaid One Fetal Doppler

High quality Fetal Doppler with fixed 2Mhz probe. Visual FHR display, excellent audio due to wide beam coverage. Made by top brand, Huntleigh.



### G2002 Fetal Doppler

Fetal Doppler with fixed 2Mhz probe. High sensitivity for audio FHR but no visual display. Detects FHR from 10th week of pregnancy.



### Huntleigh Mini Doppler

Fetal or Vascular doppler. Comes without probe which must be purchased seperately. Can choose from fetal or vascular probes. Audio only, no visual. Very clear audio quality and made by top brand, Huntleigh.



### Huntleigh Rechargeable Obstetric Doppler

Fetal doppler with fixed 2Mhz probe. High sensitivity. Audio and visual FHR. Made by top brand Huntleigh. Detects FHR from 8 weeks right through to labour. Rechargeable so no running costs.