



Welch Allyn FlexiPort™ Disposable Blood Pressure Cuffs

True Standardization in a Snap.

WelchAllyn®

Advancing Frontline Care™

Simple to order. Simple to use. Simply amazing.

Welch Allyn FlexiPort™ Disposable Blood Pressure Cuffs

True standardization in a snap.

Convenient, comfortable, and affordable, Welch Allyn disposable blood pressure cuffs are the top choice for infection-control-conscious hospitals and clinics around the world. That's because we offer virtually every type of cuff an institution may need, along with the ability to simplify inventory like no other cuff provider with our patented FlexiPort connection.

With FlexiPort Disposable Cuffs, it's easy to move patients from room to room, floor to floor, and department to department. Simply equip all devices with FlexiPort fittings, assign one cuff to your patients as they enter the hospital, and connect the patient's cuff to the device.

When you standardize on Welch Allyn FlexiPort Disposable Blood Pressure Cuffs:

- Every Cuff can be used as one- or two-tube cuff
- Part numbers are reduced by up to 60 percent, simplifying ordering and inventory
- You get every cuff size needed—from the tiniest neonatal through thigh cuffs
- Choose from either soft or vinyl fabric to fit anyone's preference
- You reduce the risk of cross-contamination

*Neonate cuffs are not FlexiPort compatible



Color coded for fast, accurate
selection of the correct cuff size



What's FlexiPort?

The FlexiPort fitting provides a single point of connection between devices and cuffs, making any cuff a one- or two-tube cuff.



Patented FlexiPort™ technology allows one cuff to work in virtually any patient care area.



Meets all the latest clinical guidelines for proper fit from AAMI and AHA

No sewn seams to tear or fray
Soft, latex-free material minimizes risk of allergic reactions and ensures patient comfort

Unique ability to connect and disconnect tubing making one- and two-tube cuffs interchangeable

Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability

REDUCE PART NUMBERS BY 60%

ONE CONNECTION POINT



The Only True Standardization Solution

With FlexiPort, Welch Allyn is the only company to offer a true cuff standardization solution for your facility. By connecting tubes at the cuff, disposable FlexiPort cuffs not only make it easier for caregivers to find and change the right cuff, but also reduce part numbers by up to 60 percent, simplifying inventory and ordering.

Contact your local Welch Allyn sales representative for more information on simplifying your cuff inventory with Welch Allyn.

Truly standardize your facility and stop worrying about tubes and connectors

- 1) Equip all blood pressure devices with FlexiPort fittings
- 2) Assign one cuff to patients upon admittance to the hospital
- 3) Have the unique ability to attach patients' cuffs to blood pressure devices—room to room, floor to floor, and department to department

The Disposable Difference

Aside from the obvious benefit of simplified ordering and application, choosing Welch Allyn FlexiPort Disposable Blood Pressure Cuffs over reusable cuffs may significantly decrease the risk of cross-contamination. Hospital-induced (nosocomial) infections strike seven to ten percent of all patients, and in addition to the negative effects on length of stay and cost, these infections are responsible for 80,000 deaths in the U.S. each year. Poorly cleaned reusable blood pressure cuffs, along with other factors, are a primary source for such bacteria.

By converting to disposable cuffs, you can significantly reduce the risk of cross-contamination and related infections.

Consider this:


Reusable blood pressure cuffs can carry such bacteria as *Klebsiella pneumoniae*, Norwegian Scabies, *E. cloacae*, and *Staphylococcus aureus*—bacteria potentially deadly to infants and adults alike.

Reused blood pressure cuffs are a primary carrier of MRSA, which is responsible for 12% of all hospital-related infections in the U.S.











Source: LONGITUDINAL EVALUATION OF NEONATAL NOSOCOMIAL INFECTIONS: ASSOCIATION OF INFECTION WITH BLOOD PRESSURE CUFF; Martin G. Myers, M.D.; (Pediatrics 1978; Vol. 61 42-45).

Source: METHICILLIN-RESISTANT STAPHYLOCOCCUS AUREUS (MRSA) AND VANCOMYCIN-RESISTANT ENTEROCOCCI (VRE) IDENTIFIED ON ENVIRONMENTAL SURFACES; J. M. Boyce, M. M. Jackson, G. Pugliese, D. B. Murray,



Simplify inventory by reducing part numbers

Once these FlexiPort fittings are attached to your blood pressure devices, you can buy cuffs without tubes and connectors—reducing part numbers by up to 60%.

		Option 1	Option 2							
		Attach above fittings directly to your current device tubing.	Buy FlexiPort cuffs with tubing and connectors to fit your blood pressure devices without making changes.							
		No Tubes or Connectors Required	One-Tube				Two-Tube			
										
Size (20 cuffs/case)	Range	FlexiPort	Screw	Tri-Purpose	Locking	Bayonet	Screw	Locking	Tri-Purpose & Empty Tube	Inflation Sys.
Small Infant (6)	7-10 cm	SOFT-06	SOFT-06-1SC	SOFT-06-1TP	SOFT-06-1MQ	SOFT-06-1HP	n/a	n/a	n/a	n/a
Infant (7)	9-13 cm	SOFT-07	SOFT-07-1SC	SOFT-07-1TP	SOFT-07-1MQ	SOFT-07-1HP	SOFT-07-2SC	SOFT-07-2MQ	SOFT-07-2TP	SOFT-07-2BV
Small Child (8)	12-16 cm	SOFT-08	SOFT-08-1SC	SOFT-08-1TP	SOFT-08-1MQ	SOFT-08-1HP	SOFT-08-2SC	SOFT-08-2MQ	SOFT-08-2TP	SOFT-08-2BV
Child (9)	15-21 cm	SOFT-09	SOFT-09-1SC	SOFT-09-1TP	SOFT-09-1MQ	SOFT-09-1HP	SOFT-09-2SC	SOFT-09-2MQ	SOFT-09-2TP	SOFT-09-2BV
Small Adult (10)	20-26 cm	SOFT-10	SOFT-10-1SC	SOFT-10-1TP	SOFT-10-1MQ	SOFT-10-1HP	SOFT-10-2SC	SOFT-10-2MQ	SOFT-10-2TP	SOFT-10-2BV
Adult (11)	25-34 cm	SOFT-11	SOFT-11-1SC	SOFT-11-1TP	SOFT-11-1MQ	SOFT-11-1HP	SOFT-11-2SC	SOFT-11-2MQ	SOFT-11-2TP	SOFT-11-2BV
Adult Long (11L)	25-34 cm	SOFT-11L	SOFT-11L-1SC	SOFT-11L-1TP	SOFT-11L-1MQ	SOFT-11L-1HP	SOFT-11L-2SC	SOFT-11L-2MQ	SOFT-11L-2TP	SOFT-11L-2BV
Large Adult (12)	32-43 cm	SOFT-12	SOFT-12-1SC	SOFT-12-1TP	SOFT-12-1MQ	SOFT-12-1HP	SOFT-12-2SC	SOFT-12-2MQ	SOFT-12-2TP	SOFT-12-2BV
Large Adult Long (12L)	32-43 cm	SOFT-12L	SOFT-12L-1SC	SOFT-12L-1TP	SOFT-12L-1MQ	SOFT-12L-1HP	SOFT-12L-2SC	SOFT-12L-2MQ	SOFT-12L-2TP	SOFT-12L-2BV
Thigh (13)	40-55 cm	SOFT-13	SOFT-13-1SC	SOFT-13-1TP	SOFT-13-1MQ	SOFT-13-1HP	SOFT-13-2SC	SOFT-13-2MQ	SOFT-13-2TP	SOFT-13-2BV

For vinyl cuffs, replace "SOFT" with "VINYL" in part number (Vinyl "TP" version not available).

Neonatal Soft Disposable Cuffs with Male Luer Slips

Size	Range	One-Tube (10 pack)	Two-Tube (10 pack)
NEO 1	3.3-5.6 cm	5082-101-1	5082-101-2
NEO 2	4.2-7.1 cm	5082-102-1	5082-102-2
NEO 3	5.4-9.1 cm	5082-103-1	5082-103-2
NEO 4	6.9-11.7 cm	5082-104-1	5082-104-2
NEO 5	8.9-15.0 cm	5082-105-1	5082-105-2
Multi-Pack	1 of each	5082-241-9	5082-241-10



For more information, contact your local Welch Allyn representative



Advancing Frontline Care™

Tel: +353 (0) 1 835 2411
Fax: +353 (0) 1 969 5009

Medguard Professional Healthcare

Email: service@medguard.ie
Website: www.medguard.ie