



## Sekusept® MultiEnzyme

# ENZYMATIC DETERGENT FOR MEDICAL DEVICES



### FAST



Sekusept Multienzyme offers excellent cleaning efficacy to help maximise the manual cleaning process and avoid time spent re-cleaning.

### SIMPLE

Our range of dispensing solutions and user friendly packaging help you to achieve the right dilution every time with minimum effort.

### SAFE

The neutral pH means Sekusept Multienzyme is both safe for the user and compatible with all medical devices. The transparent blue colour of the solution indicates that product has been added and all instruments are visible.



## Why use Sekusept Multienzyme?

- Excellent cleaning performance
- Contains a unique blend of surfactants and four enzymes to help break down organic matter
- CE Marked for use on medical devices
- Low-foaming formulation for easy rinsing and high visibility of instrumentation
- Accurate dosing with a choice of measuring beaker, pump dispenser or wall mounted dispenser
- Transparent blue colour to indicate that product has been added
- Pleasant fragrance

## Areas of use

For manual pre-cleaning of surgical instruments and medical devices including flexible endoscopes.

Also suitable for use in ultrasonic baths and automated endoscope reprocessors.

## Ingredients

Proteases, Amylases, Lipases, Surfactants, Corrosion Inhibitors, Solvents, Auxiliary agents, Emulsifiers



## Ordering information

Description	Unit size	Outer size	Ecolab re-order code
Sekusept Multienzyme P	6L	1	P304667
Sekusept Multienzyme P	2L	4	P304665
20ml pump dispenser for use with 6L bottle	20ml pump	1	TBC
10ml pump dispenser for use with 2L bottle	10ml pump	4	700064

