

Overview:

This new bariatric tilt table is designed to safely accommodate heavier patients with ease. Lowering to just wheelchair height it features electric height and tilt adjustment. The large braked castors, adjustable bariatric harnesses and patient grab handles are all designed specifically for obese patients. The tilt table is operated by a four button hand control permitting simultaneous control and fine adjustment of both height and tilt functions, with an angle indicator to display the angle.



Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Electric height & tilt operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 285kg
- 125mm heavy duty braking castors
- Hoist clearance 130mm
- Swing back patient grab handles
- Bariatric tilting divided footboard
- Battery back-up in event of power failure
- Bariatric patient harnesses
- Angle indicator

Overall Dimensions

Length:	185cm (73") approx
Width:	80cm (31") approx
Height:	48cm - 97cm approx (19" - 38") approx

Upholstery Dimensions



Recommended Options:



Patient Table




Tilting Patient Table



D Cushion

Operational Parameters



Model BARI1 - Electric version

For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Fax: +353 (0) 1 969 5009

Email: service@medguard.ie

Website: www.medguard.ie