



MEDGUARD[®]
Professional Healthcare Supplies



Chair weigher

7802

- ▶ Robust construction
- ▶ Simple weighing of persons with a mobility handicap
- ▶ Ergonomically shaped and easy to clean seat
- ▶ Hold and tare function



Features

- ▶ Ergonomically shaped and easy to clean seat
- ▶ Hinged arm- and footrests
- ▶ Energy-saving auto-off function
- ▶ Rechargeable battery operation with integrated charger
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value will be held for a few seconds after unloading
- ▶ Version .002: two rear swivel castors with rear fixing brakes
- ▶ Version .402: four swivel castors with rear fixing brakes
- ▶ Switchable to kg/lb units

Products

Ordernumber Variation	Approvable	Division g	Maximum load kg	Model
7802.70.002 Chair weigher with display Basic	No	100	200	2 rear swivel castors with rear fixing brakes
7802.70.402 Chair weigher with display Basic	No	100	200	4 swivel castors with rear fixing brakes