











A Trusted Leader in Auscultation




3M™ Littmann® Stethoscopes Product Comparison Chart

	Electronic Stethoscopes		Cardiology Stethoscopes			Classic Stethoscopes			Select/Lightweight Stethoscopes	
	Electronic Stethoscope Model 3200	Electronic Stethoscope Model 3100	Master Cardiology™ Stethoscope	Cardiology S.T.C. (Soft Touch Chestpiece) Stethoscope	Cardiology III™ Stethoscope	Master Classic II™ Stethoscope	Classic II SE Stethoscope	Classic II Pediatric and Infant Stethoscope	Select Stethoscope	Lightweight II S.E. Stethoscope
										
Overall Performance (1-10 Scale)	10+	10+	10	9	9	8	7	7	6	6
Main Application Areas	Cardiology/ High Performance Bluetooth®	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Blood Pressure, Limited Physical Assessment	Blood Pressure, Limited Physical Assessment
Chestpiece Design and Technology	Durable Metal	Durable Metal	Polished Stainless Steel Single-sided	Matte Finished Stainless Steel, Single-sided with Soft Touch Insert	Machined Stainless Steel Two-sided Tunable	Plated and Polished Alloy Single-sided	Machined Stainless Steel Two-sided	Machined Stainless Steel Two-sided	Coated Alloy Single-sided	Metal/Resin Composite Two-sided
Headset Material	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy
Standard Length	27"/69 cm	27"/69 cm	22"/56 cm 27"/69 cm	27"/69 cm	22"/56 cm 27"/69 cm	27"/69 cm	28"/71 cm	28"/71 cm	28"/71 cm	28"/71 cm
Total Weight	185 grams	185 grams	175, 185 grams	178 grams	165, 175 grams	155 grams	125 grams	105, 90 grams	115 grams	118 grams
Diaphragm Details, Diameter	51.0 mm	51.0 mm	Tunable 1.75"/4,4 cm	Tunable, Teardrop 2.15" x 1.9" 5,4 cm x 4,8 cm	Two Tunable Diaphragms 1.7"/4,3 cm	Tunable 1.75"/4,4 cm	Tunable 1.75"/4,4 cm	1.31"/3,3 cm 0.75"/1,9 cm	Tunable 1.75"/4,4 cm	Tunable, Teardrop 2.15" x 1.9" 5,4 cm x 4,8 cm
Warranty	2 Years	2 Years	7 Years	6 Years	5 Years	3 Years	3 Years	3 Years	2 Years	2 Years





3M™ Littmann® Ordering Information and Tube Color Guide



	Black	Burgundy	Navy Blue	Hunter Green	Plum	Chocolate	Gray	Pine Green	Ceil Blue	Purple	Raspberry	Caribbean Blue	Orange	Seafoam Green	Red	Lilac	Royal Blue	Ocean Blue	Peach	Pearl Pink
Model 3200	27 in. (69 cm) 3200BK27	3200BU	3200NB																	
Model 3100	27 in. (69 cm) 3100BK27	3100BU	3100NB	3100HG	3100PL															
Master Cardiology™	22 in. (56 cm) 2159																			
	27 in. (69 cm) 2160 2161 ¹ NEW! 2175² NEW! 2176³	2163	2164	2165	2167	2169	2168				NEW! 2178									
Cardiology STC	27 in. (69 cm) 4471	4472	4473		4475															
Cardiology III™	22 in. (56 cm) 3127																			
	27 in. (69 cm) 3128 3131BE ¹ 3128BRS ²	3129	3130	3134	3135	3137 3137CPR ³	3136					3138			3140					
Master Classic II™	27 in. (69 cm) 2144L 2141 ¹ 2142G ⁴	2146	2147	2632				2634	2633	2143		2630								
Classic II S.E.	28 in. (71 cm) 2201 2218BE ¹ 2201BRS ²	2211	2205	2208		2820CPR ³	2203	2818	2813	2209	2210	2206 NEW! 2823⁶	2812	2814				2819	2822	2817
Classic II Pediatric	28 in. (71 cm) 2113		2123					2131			2122	2119 NEW! 2153⁶			2113R					
Classic II Infant	28 in. (71 cm) 2114							2132			2125 NEW! 2157⁶	2124			2114R	2126				
Select	28 in. (71 cm) 2290	2293						2305	2301	2294	2296	2291				2303	2298	2306	2310	
Lightweight II S.E.	28 in. (71 cm) 2450	2451							2454			2452		2455		2453				



Medguard Professional Healthcare Supplies

Telephone: + 353 (0) 835 2411

Email: sales@medguard.ie

Fax: + 353 (0) 1969 5009

Website: www.medguard.ie

3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Master Classic II and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Bluetooth® is a registered trademark of Bluetooth SIG. Please recycle. Printed in U.S.A. © 3M 2007, 2008, 2009, 2010. All rights reserved. 70-2008-9603-6



- 1 Black Edition – black finish chestpiece and eartubes.
- 2 Brass Edition – brass finish chestpiece and eartubes.
- 3 Copper Edition – copper finish chestpiece and eartubes.
- 4 Gold Edition – gold finish chestpiece and eartubes.
- 5 Smoke Edition – smoke finish chestpiece and eartubes.
- 6 Rainbow Edition – rainbow finish chestpiece and sunshine gold eartubes.