

Ramblegard Wireless Monitor



Picture for illustrative purposes only

Product Code:

RG2

Description:

The Ramblegard wireless monitor works in conjunction with bed and chair sensor pads, floor sensor mats, or PIR, and alerts accordingly.

Features:

- Works with two pads simultaneously
- Reset button on monitor face
- Tamper proof reset option using magnetic caregiver key
- Choice of tones and volume control for alarm
- Port to connect with nurse call system (NC/NO)
- Operates on batteries or via AC adaptor
- Low battery and nurse call plug removal alert
- Pad/signal lost alert

Works in conjunction with:

- RG9 - Transmitter
- RG13 - Wireless Bed Sensor Pad
- RG11 - Wireless Chair Sensor Pad
- RG16 - Wireless Floor Sensor Mat
- RG23 - Carer Key
- RG20 - Infra Red Sensor
- RG21 - Wireless Nursecall Button

Technical Specifications:

Operating frequency: 433.92 Mhz
Operating range: 12 m
Power Supply: 6V DC, 300mA
: or 3 x AA-type alkaline batteries
Battery life: 6 months (Nominal)
Dimensions(LxWxD): 120 x 77 x 34mm
Housing material: ABS (UL-94HB, flammability grade)
Net Weight: 270g (with softcover)
Operating temp.: 5°C ~ 40°C
Rel. Humidity: < 85%

Approvals:

CE Approved & compliant with R&TTE directive
EMC: EN301489-1 & 3
Radio: EN300220-1 & 2
Alarm: EN50130-4
Safety: EN60950-1
Medical Device Immunity: EN60601-1-2
FCC Approved
FCC DOC

**It stands out from its competitors,
it is effective and it's economical.**

MEDGUARD[®]
Professional Healthcare Supplies

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Fax: +353 (0) 1 969 5009

Email: service@medguard.ie

Website: www.medguard.ie