

# Wireless<sup>📶</sup> Chair Sensor Pad

## Product Code:

RG11

## Description:

This bed sensor pad works with Ramblegard Wireless Monitors, and will trigger an alarm when someone releases pressure by exiting the chair



Picture for illustrative purposes only

## Features:

- Triggers alarm when pressure is released
- Durable and long-lasting
- Care instructions printed on pad
- Versatile and economical
- 1 year guarantee

## Works in conjunction with:

- RG2 - Ramblegard Monitor
- RG9 Transmitter

## Technical Specifications:

Cable: for 433MHz Cordless transmitter

Material: PVC & PET

Dimensions (LxWxD): 38 x 25cm

Net Weight: 100g

Operating temp.: 5°C to 40°C

Operating weight: 10 - 200 kg

Recommended Cleaning: 1:100 solution of household bleach (5.25-6.15% sodium hypo-chlorite)

**It stands out from its competitors,  
it is effective and it's economical.**

**MEDGUARD**<sup>®</sup>  
Professional Healthcare Supplies

Medguard Professional Healthcare

Tel: +353 (0) 1835 2411

Fax: +353 (0) 1969 5009

Email: [service@medguard.ie](mailto:service@medguard.ie)

Website: [www.medguard.ie](http://www.medguard.ie)