



Woburn Community 1200 Profiling Bed

Overview

Product Name: Woburn Community 1200 Profiling Bed

Product Code: HLBC791/12

Description: The Woburn Community 1200 is Harvest's first wide profiling bed, built to the same high quality as other beds in the range. Although not a bariatric bed, it has a SWL of 34 Stone. The generous width allows patients to move freely. The lockable handset allows for quick adjustments of height, back rest and knee break and the attractive vinyl-clad head and foot boards facilitate infection control.

Safe Working Load: 220kg / 34 Stone

Max. Patient Weight: 185kg / 29 Stone

Overall Dimensions: 2120 x 1320mm

Mattress Base: 2000 x 1200mm

Max. Mattress Depth: 205mm

Adjustable Height: 385 - 815mm

Angle of Back Rest: Up to 70 degrees

Angle of Knee Break: Up to 30 degrees

Bed Weight: 91.5kg (without headboards)

Key Information: 4 Section Profiling Bed
Electronically Powered Back Rest & Knee Break
Manual Leg Rest Adjustment
Vinyl Coated Head and Foot Boards
Integral Side Rails
Trendelenburg Function - Optional
Lockable Handset
Low Voltage – 29v; 9v Battery Backup
Comprehensive Manual; CE marked
Warranty: 5Yrs

Medguard Professional Health Care

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co Meath

Tel: + 353 (0) 1 835 2411

Email: service@medguard.ie

Fax: +353 (0) 1 969 5009

Website: www.medguard.ie