



Product description

Order number 75.9922.745
Product description Tube 70 ml, 55 x 44 mm, transparent polypropylene, assembled yellow cap, with standard label

Product characteristics

Label/ Print With plastic label
Graduation No
Type of cap Screw cap
Design Assembled cap, with standard label
Application Collection and storage of urine
Example of use Collection tube for urine samples

Size

Sample volume 70 ml
Max. volume 70 ml
Diameter 44 mm
Diameter of opening 44 mm
Height of product 55 mm
Length excluding cap 55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	Yellow

Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	batched

Packaging

Minimum order qty.	500
Pcs. / inner box	250
Pcs. / box	500
Pcs. / pallet	8000
Case volume	0,0889 cbm
Weight of case	7 KG
EAN of inner box	4038917152761
EAN of outer box	4038917027809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.