

**NEXT DAY
DELIVERY**

MEDGUARD®

Professional Healthcare Supplies



Rental - 24 Hour BP



ErgoScan
24 Hour ABPM System

For monitoring and diagnosis
of hypertension

Full details inside



€39.00
ex vat per month

Rental - ECG



KES-601T ECG machine

Cardiovascular diagnostics

High quality, German made

Full details inside



€69.00
ex vat per month

Rental - Vital Signs



Adviview Vital Signs Monitor

For blood pressure, temperature
& O2 saturation testing

Available in wall, desk or moveable
mobile stand models



€79.00
ex vat per month

ErgoScan 24 HR ABPM

The Ergoscan long-term blood pressure monitor stands out not only with its ease of use but also thanks to its light weight.

The long-term blood pressure monitor comes delivered with the corresponding software and a connection cable, so it can be used immediately.

Simple to install and to use, allowing easy integration of readings with patient records.

The Ergoscan long-term blood pressure monitor can be programmed to suit your surgery's individual needs.



ErgoScan 24 Hour BP Key Features:

- **PC-Based Evaluation**
Full analysis & secure data transfer
- **Attaches To Patient Files**
Readings easily attached to patient files
- **Large Internal Memory**
Memory for up to 200 measurements
- **Light & Comfortable**
Very small, only weighing 190 grams
- **Patient Friendly**
Will not over tighten on the arm
- **Configurable**
Activation of additional measurements possible

What's included?

- Ergoscan 24 HR BP Recorder
- Ergoscan ABPM Cuff - standard size (24-32cm/13cm wide)
- Ergoscan ABPM Cuff - large (32-42cm/13cm wide)
- Carrying bag
- Belt
- Charging station
- 4 NiMH rechargeable batteries
- Adult Fabricel Protective Cuff Barrier (50)
- Connecting cable to PC (USB)
- Windows software on USB

€39.00 ex vat per month
for the ErgoScan without software

€45.00 ex vat per month
for the ErgoScan with software

Kalamed KES-601T ECG Machine

The KES-601T is a small and lightweight ECG unit with 6-channel Printer. Its large LCD Touch-Display offers ease of use for checking in patient data as well as clear, accurate viewing of ECG waveforms.

The KES-601T can be used as a stand-alone ECG or connected to a PC for data transfer. Data can also be transferred via USB memory stick.

Readings can be printed directly from the machine or to an external laser printer.

ECG interpretation comes as standard, as well as an internal memory of 250 readings:



Kalamed KES-601T Key Features:

- **Large LCD Display**
6" ECG Touch-Display for 12 leads
- **Interpretative**
ECG interpretation for every reading
- **12-Lead**
12 Lead simultaneous ECG Acquisition
- **Automatic ECG Measurement**
Pattern recognition to identify irregularities
- **European Quality Standard**
Manufactured in Germany to the highest standard
- **Connective**
Readings can be transferred to patient files

What's included?

- KES-601T ECG machine
- Patient Cable with 4mm Banana Plugs - 10 Lead Diagnosis ECG Cable x 1
- Crocodile clips for ECG electrodes x 10
- Power Cord x 1
- Paper Roller x 2
- Thermo Sensitive Paper (112 mm wide x 30 m, Printing area:104 mm) x 1
- Ground Lead x 1
- Fuse x 2
- Small flat-head Screwdriver x 1
- Dust Cover x 1

€69⁰⁰ ex vat per month

Adview Vital Signs Monitor

ADC, a world leader in core diagnostic products for over a quarter of a century, have now created Adview, the first truly modular diagnostic station that grows with your needs.

Start with basic automated BP, then add thermometry and pulse oximetry as you need it.

Choose from wall mount, table mount or mobile configurations. Or convert the unit to the platform needed as your requirements change.



Adview Key Features:

- Automated Blood Pressure
Systolic, diastolic, pulse rate & MAP in 30 seconds
- “Fastemp” Technology
Oral or axillary temperature in 10 seconds
- Nellcor SP02
Accurate O2 saturation readings in 10 seconds
- Optional Integration
Electronic medical records storage
- Modular
Purchase whole or upgrade at any time
- Optional Platforms
Wall, table and mobile stands available

Choose Your Platform:

Mobile

With the five-legged mobile stand you can move the instrument from room to room for optimum workflow

Wall mount

With custom mounting plate, cuff basket included

Table

Innovative caddy that holds the unit, power supply, cuff



€79.00
ex vat per month