



# Baby scale Conformity assigned "Exklusiv"

7752

- ▶ Conformity assigned to class III, MPG class 1
- ▶ Very fine division 2 g / 5 g
- ▶ Tara function: Covers / paddings can be tared
- ▶ Hold function: The weight value will freeze after stand still



## Features

---

- ▶ Attractive baby scale with very fine division of 2g / 5g
- ▶ 5- button operation with programmable function key
- ▶ Protection Class IP X2
- ▶ Backlit LCD Display, 16 mm digit height
- ▶ Manual input of tare values by function key
- ▶ Tara function: Covers / paddings can be tared
- ▶ Hold function: The weight value will freeze after stand still
- ▶ Provisional tare
- ▶ Simple differential weighing after breastfeeding / feeding possible (BMIF)
- ▶ Analogue and digital measurements with baby height rod, optionally available
- ▶ RS 232 interface (decoupled) optional
- ▶ Units in kg
- ▶ Mains operation


## Products

---

| Order number<br>Variation            | Approvable | Division<br>g | Maximum load<br>kg | Platform size<br>mm |
|--------------------------------------|------------|---------------|--------------------|---------------------|
| 7752.01.001<br>Baby scale "Exklusiv" | Yes        | 2   5         | 6   15             | 601 x 385 x 129     |

---

## Accessories

| Order number<br>Variation |   | Description   |
|---------------------------|---|---|
| 2563.41.004               |   | Interface RS232, installed (must be included in the purchase order, no retrofit possible) |
| 5005.11.001               |  | Baby height rod   |
| 2792.20.002               |   | Alibi memory, installed (must be included in the purchase order, no retrofit possible)    |
| 2550.03.002               |   | Data cable PC/IT (Cable length 2 m)   |
| 2550.00.008               |   | Data cable PC/IT (Cable length 2 – 10 m)  |
| 2557.25.002               |   | UK adapter  |