



# Baby height rod

5005

- ▶ Digital and analogue readings
- ▶ Easy to use, quick and precise
- ▶ Hold function: The reading remains on the display
- ▶ Energy-saving auto-off function

## Features

---

- ▶ Easy to use, quick and precise
- ▶ Hold function: The reading remains on the display
- ▶ LCD display, 16 mm digit height
- ▶ Incl. adapter for attachment to baby scale 7725 or 8352
- ▶ Can also be used without a scale
- ▶ Units in cm or ft.
- ▶ Digits readable in mm,ft./inch
- ▶ Energy-saving auto-off function
- ▶ Battery operation with 2 x 1.5 V AAA batteries (included in scope of delivery)

## Products

---

Order number Variation	Approvable	Division mm	Measuring range cm
5005.11.001 Baby height rod for 7752 and 8352	No	1	35 - 81

---