

Corded Seat Mat



Product Code:

RG10

Description:

This chair sensor pad works with all Ramblegard monitors, and will trigger an alarm when someone releases pressure by exiting the chair.

Picture for illustrative purposes only

Features:

- Triggers alarm when pressure is released
- USB break-cable to protect monitor
- Durable and long-lasting
- Care instructions printed on pad
- Versatile and economical
- 1 year guarantee

Works in conjunction with:

- RG1 - Corded Ramblegard Monitor
- RG2 - Ramblegard Monitor

Technical Specifications:

Cable: 50cm coil cord with RJ10 plug

Material: PVC & PET

Dimensions (LxWxD): 38 x 25cm

Net Weight: 130g

Operating temp.: 5°C to 40°C

Operating weight: 10 - 200 kg

Internal Resistance: <1200Ω

Recommended Cleaning: 1:100 solution of household bleach (5.25-6.15% sodium hypo-chlorite)

**It stands out from its competitors,
it is effective and it's economical.**

Ramblegard[®]
stay a step ahead