


Corded Bed Sensor Mat only

Product Code:

RG12

Description:

This bed sensor pad works with all Ramblegard monitors, and will trigger an alarm when someone releases pressure by sitting up or exiting the bed



Picture for illustrative purposes only

Features:

- Triggers alarm when pressure is released
- USB break-cable to protect monitor
- Durable and long-lasting
- Care instructions printed on pad
- Versatile and economical
- 1 year guarantee

Works in conjunction with:

- RG1 - Corded Ramblegard Monitor
- RG2 - Ramblegard Monitor

Technical Specifications:

Cable: 2.5m breakaway cord with RJ10 plug (N/O)

Material: PVC & PET

Dimensions (LxW): 76 x 51cm

Net Weight: 380g

Operating temp.: 5°C to 40°C

Operating weight: 10 - 200 kg

Internal Resistance:<1200Ω

Recommended Cleaning: 1:100 solution of household bleach (5.25-6.15% sodium hypo-chlorite)

**It stands out from its competitors,
it is effective and it's economical.**

Ramblegard
stay a step ahead