

Arm and leg rest on castors
Screw height adjustment

S2112027

+



Product description

Material:	steel, chromium plated
Material of foam:	PU foam
Regulating range:	appr. 680 - 880 mm
Foam pad:	appr. 300 x 100 mm
Base diameter:	appr. 635 mm
Castors:	twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable
Total weight (net):	appr. 3,0 kg



Arm and leg rest on castors

Screw height adjustment



Application

- + optimal hygienic characteristics (wipe disinfection)
-
- + very easy and fast assembly through a plug in and out system with only 1 screw connection

Technology

- + castors with non-marking tread

Economy

- + complete recyclability