

Arm and leg rest on anti-slide stoppers
Screw safety height adjustment

S2102037

+

made
in
Germany

Product description

Material:	steel, chromium plated
Material of foam:	PU foam
Regulating range:	680 - 880 mm
Foam pad:	300 x 100 mm
Base diameter:	appr. 635 mm
Total weight (net):	appr. 3,0 kg



Arm and leg rest on anti-slide stoppers

Screw safety height adjustment



Application

- + optimal hygienic characteristics (wipe disinfection)
-
- + very easy and fast assembly through a plug in and out system with only 1 screw connection

Economy

- + complete recyclability