

Arm and leg rest on castors  
single handed safety height adjustment

S1112017

+



Product description

Material:	steel, chromium plated
Material of foam:	PU foam
Regulating range:	680 - 880 mm
Foam pad:	300 x 100 mm
Base diameter:	appr. 635 mm
Castors:	twin castors Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable
Total weight (net):	appr. 3,0 kg



# Arm and leg rest on castors

## single handed safety height adjustment

---



---

### Application

- + optimal hygienic characteristics (wipe disinfection)
- 
- + very easy and fast assembly through a plug in and out system with only 1 screw connection

### Technology

- + castors with non-marking tread

### Economy

- + complete recyclability