



MEDGUARD®

Professional Healthcare Supplies

	Susie Simon® Nursing Care Simulator	Susie Simon® Nursing Care Simulator	Simple Simon® Nursing Care Simulator	Simple Simon® Nursing Care Simulator	3B Scientific® Patient Care Manikin PRO
Item No.	without ostomy: 1017535	with ostomy: 1005785	without ostomy: 1005803	without ostomy: 1005807	1018816
Gender	male and female	male and female	female	male	male and female
Wig for hair care	✓	✓	✓	-	✓
Movable jaw with tongue	✓	✓	✓	✓	✓
Removable dentures	✓	✓	✓	✓	✓
Ear irrigation	✓	✓	✓	✓	✓
Nasogastric lavage and gavage	✓	✓	-	-	✓
Tracheotomy care	-	-	-	-	✓
IM sites at upper arm and in buttock	✓	✓	✓	✓	✓
IM site at the upper legs	-	-	-	-	✓
Enema administration	✓	✓	-	-	✓
Catheterization	✓	✓	-	-	✓
Realistic vagina and cervix	✓	✓	✓	-	-
Vaginal douching and pap smear	✓	✓	-	-	-
Sculpted stomas	-	✓	✓	-	✓
Surgical & amputation wound	-	-	-	-	✓
Diabetic foot & decubitus ulcers	-	-	-	-	✓
Laceration wound	-	-	-	-	✓
Detachable and removable internal tanks	✓	✓	-	✓	✓

HPRA Wholesalers authorisation No. W00426/00001

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Unit 28B Ashbourne Business Centre, Ashbourne, Co Meath, A84 WA49

Transverse colostomy, ileostomy, and suprapubic stoma	-	✓	-	-	only colostomy
---	---	---	---	---	----------------