



Laundry Bags


Safetex Bags

Safetex is a re-usable, self opening laundry bag, designed for the collection of heavily soiled and infested linen. Safetex has unique welded seams to ensure the bag is totally fluid and odour proof. Safetex opens in the wash process once it is safely inside the washing machine. Safetex bags are environmentally friendly and a cost effective alternative to single use disposable red bags, having a life span of up to 1000 processings. Available in red (for infested linen) or blue and green (for delicate items that cannot withstand a boil wash).

Safetex Bags			
	Product Code	Description	Size
	LB5Q/SL	Snap Closure	60 x 75 cm



Procedure For Use



Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Fax: 01 969 5009

Website: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, A84 VK52