

**Sterile
Examination Gloves**

MEDGUARD®



Sterile Nitrile Exam Gloves

Sterile, latex free, Nitrile gloves - Ready to use - Powder-free
Suitable for high risk environments - Individually wrapped pairs



ORDER NOW

01 835 2411
orders@medguard.ie
www.medguard.ie

Inspiration, Innovation, Integrity

- Glove Information: Sterile, Nitrile, Powder-free, Light Blue Color and Ambidextrous
- Single chlorinated to provide exceptional dry hand donning capability
- Nitrile formulation promotes better puncture and chemical resistance properties
- Non-Latex formulation eliminates exposure to Type I latex protein sensitivity
- Finger tips/full textured with excellent working grip to enable efficient instrument handling
- Beaded cuff design to don gloves efficiently and securely
- Undetectable Chemical Residue to minimize Type IV hypersensitivity
- Suitable for all types of non-invasive sterile procedures
- High durability against snags and tears to provide the maximum protection

Tests performed:

- Sterilized by Gamma Irradiation (ISO 11137)
- Biocompatibility Tested (ISO 10993)

Ordering Information

INTW82ES	Small (50 Pairs)
INTW82EM	Medium (50 Pairs)
INTW82EL	Large (50 Pairs)
INTW82EXL	Extra Large (50 Pairs)

Outer case available - 6 boxes of 50 Pairs

Issue 202047

Unit 6/7, Block 4, Ashbourne Business Park, Co Meath t: 01 835 2411 e: orders@medguard.ie