

Specification Sheet - Medguard Latex Powder Free Examination Gloves (10 x 100)

Material	Natural High Grade Rubber Latex
Colour	Natural
Design & Feature	Ambidextrous, palm textured, beaded cuff
Powdered	No
Sterile	No
Type	Polymer Coated
Latex Free	Yes



Characteristics	Inspection Level	Acceptable Quality Level	Reference Standard
Dimensions	Median of 13 Test Pieces	Median of 13 Test Pieces	EN455 - 2 : 2015
Physical Properties	Median of 13 Test Pieces	Median of 13 Test Pieces	EN455 - 2 : 2015
Freedom From Holes (Water Tight Test)	GI	1.5	EN455 - 1 : 2000
Visual Defects: Major Visual Minor Visual	GI	2.5 4.0	In House Practice
Packaging Defects: Regulatory Visual Critical (incl. Gloves Counting)	GI GI S2	Unacceptable at any level 4.0 4.0	In House Practice
Powder Free Residue	N=5	-	EN455 - 3 : 2015 ASTM D3578 - 05 (2015) ASTM D6124 - 06 (2017)
Mix Size / Mix Glove / Mix Hand	Not Allowed	Not Allowed	Not Allowed

Description	Standard	
	Before Aging	After Aging
Median Force At Break (N)	Min 6 (Typical Value: 6 - 8)	Min 6 (Typical Value: 6 - 7)



Specification Sheet - Medguard Latex Powder Free Examination Gloves (10 x 100)

Description	Size	Standard
Length (mm)	All Sizes	Min 240
Palm Width (mm)	XS S M L XL	76+/-3 84+/-3 94+/-3 105+/-3 113+/-3
Thickness (mm) Single wall	All Sizes	Finger: 0.15+/-0.02 (Typical Value: 0.14 - 0.16) Palm: 0.12+/-0.02 (Typical Value: 0.12 - 0.14)



Key Benefits:

- Low protein levels decrease the potential for allergic reactions
- Highest quality powder-free latex gloves, stronger than the standard quality
- 6 gram latex offers a safer work environment with reduced puncture risk
- Natural in colour
- Ideal for use in medical examinations and procedures
- Made from Natural Rubber Latex
- Available in sizes Small to X-Large
- Supplied per case of 1000 (10 boxes of 100)

Product Codes

CS-6001	Extra Small (100 x 10)
CS-6002	Small (100 x 10)
CS-6003	Medium (100 x 10)
CS-6004	Large (100 x 10)
CS-6005	Extra Large (100 x 10)