



Specification Sheet - Disposable Aprons On A Roll - Blue & White

CS-16000 White Disposable Aprons On A Roll (200 x 5)

Material	LDPE	Thickness	16mu	Ignition Point	340°C
Length	107cm	Density	0.92g/cm ³	Toxicity	Non-Toxic
Width	69cm	Softening Point	90-100°C	Characteristic	Acid & alkali resistance

CS-17000 Blue Disposable Aprons On A Roll (200 x 5)

Material	LDPE	Thickness	16mu	Ignition Point	340°C
Length	107cm	Density	0.92g/cm ³	Toxicity	Non-Toxic
Width	69cm	Softening Point	90-100°C	Characteristic	Acid & alkali resistance

