



Specification Sheet - Mob Caps

32000 Mob Caps - Blue

Material	Polypropylene nonwoven	Weight	10 gsm	Colour	Blue
Sterile	No	Disposable	Yes - Single use	Size	21"

- Used in hospital, clinic, pharmacy, laboratory, food processing, beauty salon, electronics, industry, etc.
- Latex & fiberglass free, Non-irritating to skin, soft, comfortable, odorless
- Double elastic fixture for secure fitting
- Quality Systems: ISO 13485, ISO 9001 , CE
- Shelf Life: 3 years
- Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light
- Fine texture to keep the dust and bacteria away
- Non-toxic, environmental, odourless and tasteless
- Comfortable, soft and lightweight

Order Code	Description
32000-SINGLE	Mob Caps Blue (100)
CS-32000MG	Mob Caps Blue (10 x 100)

