

Creating a  
Safe Environment

# MEDGUARD®




## Covid-19 Risk Prevention

High filtration surgical masks that provide protection & comfort



High Filtration  
Capacity



Comfortable  
Nose Pad



Fluid  
Resistant



Soft Elastic  
Earloop

### Medguard Surgical Facemasks With Earloops

#### Features:

- Low resistance to breathing, light & comfortable
- Non-sterile, single-use
- Environmentally friendly
- Elastic earloop with dural point attachment helps provide a secure seal
- High filtration capacity
- Soft nose pad for easy and comfortable fit
- Adjustable nosepiece for extra comfort
- Fluid resistant
- Length 17.5cm, width 9.5cm

#### Conformity Standards/Regulation:

- Class I Medical Device as per MDD 93/42/EEC (Conform to EN 14683:2019)
- Produced under EN ISO 13485:2016 Quality Management
- Produced under ISO 9001:2015 Quality Management

#### Specifications:

- BFE (Bacterial Filtration Efficiency)  $\geq 98\%$
- Breathability (Delta P)  $< 49.0\text{Pa}/\text{cm}^2$
- Microbial Cleanliness  $\leq 30\text{cfu}/\text{g}$
- Splash resistance pressure 16kPa
- Shelf life 3 years

**ORDER NOW**

01 835 2411  
orders@medguard.ie  
www.medguard.ie

Order Number

23000R-SINGLE

Description

MG Facemasks With Loops (50)

Issue 202040

Unit 6/7, Block 4, Ashbourne Business Park, Co Meath t: 01 835 2411 e: orders@medguard.ie