

Specification Sheet - Medguard Surgical Masks

23000R Medguard Surgical Facemasks

Specifications:

- BFE (Bacterial Filtration Efficiency) $\geq 98\%$
- Breathability (Delta P) $< 49.0 \text{ Pa/cm}^2$
- Microbial Cleanliness $\leq 30 \text{ cfu/g}$
- Splash resistance pressure 16kPa
- Shelf life 3 years

Conformity Standards/Regulation:

- Class I Medical Device as per MDD 93/42/EEC (Conform to EN 14683:2019)
- Produced under EN ISO 13485:2016 Quality Management
- Produced under ISO 9001:2015 Quality Management

Features:

- Low resistance to breathing, light & comfortable
- Non-sterile, single-use
- Environmentally friendly
- Elastic earloop with dural point attachment helps provide a secure seal
- High filtration capacity
- Soft nose pad for easy and comfortable fit
- Adjustable nosepiece for extra comfort
- Fluid resistant
- Length 17.5cm, width 9.5cm

Outer Layer	Spun-bond polypropylene 25g/m ²
Filter Layer	Melt-blown filter polypropylene 25g/m ² , BFE 98%
Inner Layer	Spun-bond polypropylene 25g/m ²
Other Base Material	PVC coated wire
Colour	Blue
Size	17.5cmL x 9.5cmW

Order Number	Description
23000R-SINGLE	MG Facemasks With Loops (50)

