

## Specification Sheet - Medguard Vinyl Powder Free Examination Gloves Blue (10 x 200)

Material	Vinyl
Colour	Blue
Design	Ambidextrous
Powdered	No
Sterile	No
Latex Free	Yes



### Particulars:

Standard	Reference
Specification	Smooth Surface –PVC
Standards conform to	EN 455 Part 1, 2, 3 & 4
Conformity Performance	MDD 93/42/EC
Requirement Performance	Essential Requirements
Safety Performance	Risk Validation
Certification of Manufacture	ISO 9001: certified QA Ref.: 60251

### Sampling plan:

Defects	Inspection Level	AQL
Pinholes	Gen. Ins. Level I	1,5
Visual defects (major)	S-4	2,5
Visual defects (minor)	S-4	4,0
Physical Dimensions	S-2	4,0
Physical Properties	S-2	4,0

### Physical dimensions in mm and properties:

Size	Length	Width of palm	Weight (g)	Fingertip	Palm Thickness	Cuff
X-Small	240 ± 10	80 ± 5	3.5-3.7	0.06 ± 0.03	0.06 ± 0.03	0.06 ± 0.03
Small	240 ± 10	85 ± 5	3.8-4.0	0.06 ± 0.03	0.06 ± 0.03	0.06 ± 0.03
Medium	240 ± 10	95 ± 5	4.3-4.5	0.06 ± 0.03	0.06 ± 0.03	0.06 ± 0.03
Large	240 ± 10	105 ± 5	4.8-5.0	0.06 ± 0.03	0.06 ± 0.03	0.06 ± 0.03
X-Large	240 ± 10	110 ± 5	5.0-5.2	0.06 ± 0.03	0.06 ± 0.03	0.06 ± 0.03
Elongation						Min 300%

## Specification Sheet - Medguard Vinyl Powder Free Examination Gloves Blue (10 x 200)



Dimensions L x W x H (mm)	Packing per unit	Stacking per unit
Pinholes	Gen. Ins. Level I	1,5
Visual defects (major)	S-4	2,5
Visual defects (minor)	S-4	4,0
Physical Dimensions	S-2	4,0
Physical Properties	S-2	4,0

### Key Benefits:

- Best selling premium quality vinyl powder-free gloves
- Blue in colour
- Exceptional level of sensitivity
- Comfortable to wear with excellent fit
- Latex-free and therefore do not cause any skin allergies or reactions
- Available in sizes Small to Extra Large
- Supplied per case of 2000 (10 boxes of 200)

### Product Codes

CS-5002	Small (200 x 10)
CS-5003	Medium (200 x 10)
CS-5004	Large (200 x 10)
CS-5005	Extra Large (200 x 10)