

Material Safety Data Sheet Romed Alcohol Preps



IDENTIFICATION

Product name: Alcohol preps
Ref.: PREP-2000

PURPOSE

Antiseptic for cleaning skin before injection.

COMPOSITION

Prep: non woven material (mixture of polyester and viscose), 50 g/m²
Thickness 0,035 mm, 2-ply,
Dimensions prep: 65 x 30 mm ± 2 mm
Active ingredient: 70% isopropyl alcohol
Inactive ingredient: water
Liquid 0,4 – 0,6 ml per package
Alcohol contents ± 0.02 g/cm², which is ± 0.45 per piece

Foil pacakage: paper, aluminum foil and plastic. weight single pack: 73 g/m², thickness: 0.06 mm

USE

For preparation of the skin prior to an injection. Wipe injection site vigorously.
Store at room temperature
Do not use with electrocautery procedures.

NORMAL PRECAUTIONS

Avoid contact with eyes. If contacted, flush eyes with copious quantities of water.
If taken in, DO NOT induce vomiting, but obtain medical attention immediately.
Stop use if irritation and redness develop.
Ask a doctor if condition persists for more than 72 hours

OTHER PRECAUTIONS

Single use, discard after use. Keep at room temperature. For external use only. Flammable, keep away from fire or flame. Keep out of reach of children.

Packing: 100 pieces each packed in paper pouch in an innerbox, 50 innerboxes in an export carton.