

PRODUCT SPECIFICATION

Product name	Mediguard Apron
Range	Universal
Product Ref	UN440050
Product Description	Disposable Apron with a halter neck and waist ties
Marking	Product name, Universal Logo, Code Ref and single use symbol, CE Mark (minimal risk)
Material	Polyethylene 50 microns
Colour	Green
Dimensions	Length 150cm Width 91cm
Sterilisation	Non-Sterile



Primary use	For use in areas requiring protection against low level hazards, the effect of which are superficial and minimal
--------------------	--

Packaging	MOQ 200 pcs Aprons are individually wrapped 25 Aprons per Inner Box / 8 Inner boxes per carton
------------------	--

Conformity to Regulation	(EU) 2016/425 Minimal Risk PPE Category 1
Storage conditions	Store in a cool dry place away from heat and direct sunlight
Shelf life	5 years if stored as recommended above
Additional information	The product is latex free and non-sterile.

REVISION:	1.0	DATE	13/03/2020	AUTHOR	A Gidley
------------------	-----	-------------	------------	---------------	----------

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare T: +44 (0) 1530 278 111 W: www.365healthcare.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
---	--