



MEDGUARD®

Professional Healthcare Supplies

Rock Pin
Spinal board



The spinal board Spencer Rock is a device for recovery, immobilization, and transfer to the stretcher of patient with suspected spinal injuries.

Specific characteristics

- The number of handles, their positioning and their shape, allow a firm grip even in difficult conditions
- The tapered shape allows to save space and to place the device inside Spencer basket stretchers
- Dedicated holes and pins to fix belts and Spencer head immobilizers
- Smooth and filleted surfaces to allow easy cleaning
- Robust and lightweight structure to ensure impact resistance
- Allows preliminary X-Ray diagnostics

Technical data

Dimensions	1840x445x50 ± 10 mm
Materials	PE, PU
N° of handles	16
N° of pins	14
Loading capacity	200 kg
Weight	7,3 ± 0,1 kg

Optionals

ST02105A	10g vertical fastener for spinal boards
ST02108A	10g horizontal fastener for spinal boards
ST02015A	PIN STRAPS- SYSTEM OF BELTS WITH HOOKS
ST02035A	Rock straps - straps system for spinal boards
SH00200A	Spencer contour – head immobilizer
ST02106A	Mark Up – system for perm.personalisation

Class I Med.Device according to Dir. 93/42/CEE

ROCK PIN SPINE BOARD YELLOW W/PINS	
ST02010B	
CND Classification	V0804
Registration number	101223

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Fax: 01 969 5009

Website: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, A84 VK52