

WOBURN COMMUNITY LOW PROFILING BED



Product Code:	HLBC791L HLBC791LAR (with Auto-Regression)
Description:	<ul style="list-style-type: none">• Built to the same specification as our popular Woburn Community Bed• Offers the benefit of a lower entry and exit height• Ideal for many environments, fully compliant with legislation.• Comfort adjustments of height, back rest or knee break done simply using the lockable handset.
Safe Working Load:	220kg / 34 Stone
Max. Patient Weight:	185kg / 29 Stone
Overall Dimension:	2115 x 1015mm
Mattress Base:	2000 x 900mm
Max. Mattress Depth:	205mm
Adjustable Height:	220mm-640mm
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Break:	Up to 30 degrees
Bed Weight:	81kg + 11kg with rails
Key Information:	4 Section Profiling Bed Electronically Powered Back Rest & Thigh Rest Manual Leg Rest Adjustment Wooden Head and Foot Boards Side Rails - Optional Auto-Regression and Inbuilt Length Extension - Optional Trendelenburg Function as standard Lockable Handset Low Voltage – 29v; 9v Battery Backup Comprehensive Manual Warranty: 5yrs

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Fax: 01 969 5009

Website: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, A84 VK52