



Model 202

Treatment Couch

MEDGUARD®

Professional Healthcare Supplies

Overview:

Our great value, strong two section examination couch meets the basic needs of most clinics and private practices. A wide range of options are available to expand the use of this couch.

Standard Features

- UK Manufactured
- 2 year parts guarantee
- Electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 200kg
- Toe braked castors
- Hoist clearance 130mm
- Breathing hole
- Paper roll holder




Overall Dimensions

Length: 188cm (74") approx
 Width: 63cm (25") approx
 Height: 51cm - 94cm approx
 (20" - 37") approx

Upholstery Dimensions



Recommended Options:



Arm Tables



Foot Switch



Folding Cotside





Stirrups



Knee Troughs

Operational Parameters



Model 202E - Electric version

For Your Peace of Mind

- 2 year parts guarantee
(accidental damage to upholstery excluded)
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5
(BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982,
Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermiabile, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, cream/white will be supplied.



Cream/White



Dark Blue



Mid Blue

MEDGUARD[®]
Professional Healthcare Supplies