


Overview:

Flexible phlebotomy chair with unrivalled patient comfort and access. The easy removal of the fully adjustable arm rests mean this chair can also be used for general examination. Optional accessories can also be added to expand the use of this chair.




Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Hydraulic or electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Phlebotomy arm rests
- Hoist clearance 130mm


Overall Dimensions

Length:	177cm (70") approx
Width:	63cm (25") approx
Height:	52cm - 101cm approx (20" - 40") approx

Upholstery Dimensions



Recommended Options:



Stainless Tray



Battery Back-up



D Cushion



Motor Back




CPR Motor Back



Paper Roll Holder

Operational Parameters



Model 93BE - Electric version

Model 93BH - Hydraulic version

For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



HPRA Wholesalers authorisation No. W00426/00001

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Unit 28B Ashbourne Business Centre, Ashbourne, Co Meath, A84 WA49