

# Evac+Chair® 700H

FOR ASCENDING AND DESCENDING



ASCENDS  
AND  
DESCENDS



The Evac+Chair® 700 Transeat allows two people to safely transport a patient up and down stairs without carrying or heavy lifting and is well suited to spiral stairs and where space is restricted.

## Features

### 1 Telescopic extendable handle

For comfortable and safe operator positioning when encountering different passenger sizes

### 2 Stainless steel frame construction

Increased strength and robustness for continued use

### 3 Adjustable brake

For controlled descent whatever the passenger weight range

### 4 Detachable seat and backrest

Easy to clean and simple to replace if damaged or worn

### 5 Adjustable patient belts


Waist/chest, and leg belts for increased passenger safety

### 6 Front guide handles

For more controlled ascent and descent

### 7 Footrest

For passenger comfort during transportation



## Includes

- LIFETIME WARRANTY
- Patient Restraints (Chest & Foot)
- Wall Mounting Brackets
- Dust Cover
- Integrated handles with no pinch points
- Instructional DVD

## Benefits

- FDA Registered
- Easy to use
- Suitable for both ascent and descent
- Quick and easy to open
- Can be used on spiral staircases
- Compact storage
- No carrying required
- Reduces potential injury and claims
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act
- Wall mountable



Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	43"	18"	14"	32lbs	350lbs

VIDEO LINK: [evac-chair.com/watch-video/](http://evac-chair.com/watch-video/)

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®. Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11, EXCLUSIVE LIFETIME WARRANTY.

Medguard Professional Healthcare

Tel: 01 835 2378

Email: [orders@medguard.ie](mailto:orders@medguard.ie)

Fax: 01 969 5009

Web: [www.medguard.ie](http://www.medguard.ie)

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, Ireland