



LUXIFLEX LED

Luxiflex is a cold light luminaire, versatile, of compact design, special for the examination, with great benefits.

Flexibility

Flexible arm height adjustable 85 cm.
It offers a wide reaching, maximum comfort and flexibility.

Compact design

Luxiflex offers easy use, thanks to its compact LED head 8 cm in diameter and the switch located at the base of the lamp, makes this model an ideal luminaire for the most demanding examinations.

Mounting options

Standard version includes AH table clamp.
As accessories, we offer a wide selection of supports such as wall clamps, wall extension arm and trolley stand.

High performance LED

Luxiflex LED 15.000 lux at 50 centimeters.

LUXIFLEX LED

Technical Details (Cod. I0104AOLED)	
Light source	LED 6W 15°
Lamp lifetime	50.000 hours
Lux a 50 centimeters	15.000 lux
Colour temperature	4.500 °K
Colour rendering	80%
Arm length	72 cm
Colour	White





Wall support
Cod. B0101800



AH table clamp
Cod. B010190



Rail Plus bracket
Cod. 82026



Extension wall arm
Cod. B0101S0



Trolley stand
4 kg Cod. 09655
8,8 kg Cod. 09651

Medguard Professional Healthcare

Tel: 01 835 2378

Email: orders@medguard.ie

Fax: 01 969 5009

Web: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, Ireland