



# CHECK4<sup>®</sup> HCG (HL)

**Low sensitivity (1000 mIU/mL) hCG rapid test to  
confirm termination of pregnancy given under  
medical supervision**

---

Medguard Professional Healthcare

Tel: 01 835 2378

Email: [orders@medguard.ie](mailto:orders@medguard.ie)

Fax: 01 969 5009

Web: [www.medguard.ie](http://www.medguard.ie)

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, Ireland

## CHECK4® HCG (HL)

Product Code: 3066-1000AG

Check4® hCG (hl) is a rapid, qualitative test to check the outcome of first trimester Termination of Pregnancy (ToP) or miscarriage. In 5 to 10 minutes, this test will detect an hCG level of 1000 mIU/mL or more.



QUADRATECH DIAGNOSTICS

Check4



### TRUSTED

- Market leader chosen by renowned professionals
- Manufactured in the EU

### QUICK

- Results in just 5 to 10 minutes

### COST EFFECTIVE

- Convenient and reduces follow up time

### ACCURATE & RELIABLE

- Check4® hCG (hl) can detect hCG >1000 mIU/mL in over 99% of cases
- CE marked for home use



### HOW DOES IT WORK?

The Check4® hCG (hl) test specifically detects the hCG hormone in urine to check whether the hormone level is 1000 mIU/mL (determined against World Health Organisation reference) or greater.

If treatment for Termination of Pregnancy or miscarriage has been successful then the hCG level should drop below 1000 mIU/mL within 10 to 20 days. Check4® hCG (hl) will not detect hCG levels lower than 1000 mIU/mL, while classical pregnancy tests will be positive.

*Check4® hCG(hl) is not a classical pregnancy test and cannot be used for pregnancies where the initial hCG was lower than 1000mL.*