


PRODUCT SPECIFICATION SHEET (EN standard)

Reference code	Brand name	Product code	System code	SAP no.	Revised on
INTW80Exx	Intouch V Natural	OLPYT-ST5	PYT-240-ST5	ER36801CA-NA	1 Apr 2019
Product information					
Type	Sterile powder free latex examination gloves				
Material	Natural rubber latex				
Glove colour	Natural white				
Glove design	Ambidextrous, beaded cuff and full textured surface				
Surface treatment	Polymer coating on donning surface				
Sterilization method	Gamma irradiation				
Shelf life	5 years upon manufactured date				
Storage	Store under cool and dry conditions. Avoid direct sunlight, heat and petroleum products.				
Country of origin	Malaysia				
Physical dimension					
Glove size	Palm width		Glove length		
	EN455-2 standard	Product specification	EN455-2 standard	Product specification	
XS	≤ 80 mm	75 ± 5 mm	≥ 240 mm	≥ 240 mm	
S	80 ± 10 mm	85 ± 5 mm	≥ 240 mm	≥ 240 mm	
M	95 ± 10 mm	95 ± 5 mm	≥ 240 mm	≥ 240 mm	
L	110 ± 10 mm	106 ± 5 mm	≥ 240 mm	≥ 240 mm	
XL	≥ 110 mm	116 ± 5 mm	≥ 240 mm	≥ 240 mm	
Single wall thickness					
	EN455 standard	Product specification			
Finger thickness	-	≥ 0.14 mm			
Palm thickness	-	≥ 0.12 mm			
Cuff thickness	-	≥ 0.10 mm			
Physical properties					
	Condition	EN455-2 standard	Product specification		
Force at break	Unaged	≥ 6 N	≥ 6 N		
	Aged	≥ 6 N	≥ 6 N		
Residue content					
	EN455-3 standard	Product specification			
Protein content	≤ 50 µg/g	≤ 50 µg/g			
Powder residue	≤ 2 mg/glove	≤ 2 mg/glove			
Pre-shipment inspection					
	Inspection level	Acceptance quality limit (AQL)			
Freedom from holes; water tightness test (EN455-1)	G1	1.5			
Visual inspection	Critical defects	G1	Accept no defect		
	Major defects	G1	2.5		
	Minor defects	G1	4.0		
Physical dimension	S-2	4.0			
Physical properties	S-2	4.0			
Protein content	N = 4	N/A			
Powder residue	N = 5	N/A			
Quality certifications					
ISO 9001:2015; ISO 13485:2016; EN 13485:2016 & EN ISO 13485:2016					
Prepared by			Company stamp		
Name : Dr. Catherine Voo Designation : Technical Product Specialist Signature : 					