

MEDGUARD®

Professional Healthcare Supplies



Mitella Comfort

Stable arm support

Mitella Comfort

Stable arm support

Actimove® Mitella Comfort is a stable arm support. Patients will benefit from the easy application and breathable material.

- The 5.5cm wide sling eliminates the need for a moving neck pad, while offering comfort by reducing pressure on the back of the neck
- Soft thumb loop allows a stable and supported position of the forearm whilst still allowing full finger movement
- Simple to apply and adjust for an individual fit, O-rings allow simple application and adjustment
- One-handed, self-application is achievable, aided by the hook and loop fastenings and the open and close pocket design
- The gentle breathable pocket material provides protection for clothes or bandaging from dirt and moisture

Indications:

- Post-operative support
- Fracture management such as clavicle, scapular, humeral, and greater tuberosity fractures
- Joint dislocations such as acromioclavicular and sternoclavicular dislocations and shoulder luxation
- Inflammatory diseases such as calcific tendinitis of the rotator cuff
- Minor instability of the shoulder
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilisation in emergency or paramedic care



- Machine washable at 40°C
- Radiolucent
- Metal-free
- Latex-free



Product Code	Size	Measurements	Pack Quantity
72819-61	Small	< 37cm	1 x 1
72819-62	Medium	37cm - 40cm	1 x 1
72819-63	Large	40cm - 43cm	1 x 1
72819-64	Extra Large	> 43cm	1 x 1

Medguard Professional Healthcare



Actimove®,
an Essity brand



Tel: 01 835 2378

Email: orders@medguard.ie

Fax: 01 969 5009

Web: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, Ireland