


M95

## Simulator Pediatric IV Hand



### Set of pediatric hands for practicing IV

Become confident in safe insertion of intravenous cannulas on small infant hands with thin and soft veins

#### FEATURES

- 1 | Fingers and dorsum of hand are flexible, able to be gripped, and repositionable as in real-life procedures.
- 2 | Two puncture sites on the back of the hand
- 3 | Practice dressing and securing the cannula

#### SKILLS

- | Pediatric dorsal veins injection
- | Blood collection
- | Intravenous infusion

#### DESCRIPTIONS

##### SPECIFICATIONS

Size: W18 x D30 x H9.5 cm / W7.1 x D11.8x 3.7 inch (one-year-old)  
W18 x D36 x H10 cm / W7.1 x D14.2x 3.9inch (three-year-old)

##### SET INCLUDES

2 injection arm-hands (one-year-old)	1 spoon
2 injection arm-hands (three-year-old)	1 plastic jar
2 stands	1 plastic beaker
1 circulation pump	1 storage case
1 cotton swab (red)	1 instruction manual

Power: AC100-240V  
Power Consumption: 5W

##### MATERIALS

Soft resin, hard resin  
Latex free

##### REPLACEMENT PARTS

11350-010	1 replacement injection arm-hand (three-year-old)
11350-020	1 replacement injection arm-hand (one-year-old)
11388-400	10 cotton swab for coloring (red)

##### RECOMMENDED DEVICES

indwelling needle: 23G  
butterfly needle: 23G

\* Needles are not included.

M43D


## Pediatric Lumbar Puncture Simulator II

### Acquire skills for safe pediatric lumbar puncture

The full body is soft and bendable to facilitate realistic body positioning required for a safe and successful lumbar puncture.



Level of posterior superior iliac crests



#### Product Supervision

Takahiro Amano, M.D.  
Vice President  
Senior Vice-Dean, Postgraduate School  
Professor and Director, Center of Postgraduate  
Medical Education International University of Health  
and Welfare,  
Honorary Director, Sanno Medical Center  
Gregory A. Plotnikoff, M.D., M.T.S., F.A.C.P.  
Medical Director, Abbott Northwestern Hospital

#### FEATURES

- 1 | Soft, seamless baby manikin allows practice in safe patient positioning and holding.
- 2 | Two types of puncture pads: deep or shallow
- 3 | Realistic needle tip resistance
- 4 | Anatomically correct landmarks for palpation

#### SKILLS

- | Patient positioning
- | Palpation of landmarks
- | Lumbar puncture and CSF collection
- | CSF pressure measurement

#### DESCRIPTIONS

##### SPECIFICATIONS

Size: W43 x D34 x H20 cm  
W17 x D13.4 x H7.8 inch

Weight: 1.5 kg / 3.3 lbs

##### SET INCLUDES

1 baby manikin (7-10 months)	1 syringe (50mL)
2 pediatric LP block (shallow)	1 instruction manual
2 pediatric LP block (deep)	
1 irrigation bag	
1 stand for irrigation bag	

##### MATERIALS

Soft resin  
Latex free

##### REPLACEMENT PARTS

11348-100-04	stand for irrigation bag
11348-100-05	irrigation bag
11348-410	pediatric LP block (shallow)
11348-420	pediatric LP block (deep)

##### RECOMMENDED DEVICES

22G-24G spinal needle

\* Needles are not included.