

## Material Safety Data Sheet

### Section 1: Manufacture's Information

Manufacture's name: Blue Sail Medical Co., Ltd  
Adress: No.48 Yinuo Road,Linzi District,Zibo City,Shandong Province,China 255414  
Tel. 0086-533-7480107 Fax. 0086-533-7486199

### Section 2 : Product Information

Product Name:	DISPOSABLE FACE MASK(TYPE IIR)
CAS#	Not Established
LAT/IACO Hazard Rating	Not Restricted
D.O.T. Rating	Not Restricted

### Section 3 :Physical And Chemical Properties

External ply: SPP made from polypropylene(CAS#9003-07-0) Colour: Blue  
Internal ply: SPP made from polypropylene(CAS#9003-07-0) Colour: White  
Middle ply: Meltblown fabric made from polypropylene(CAS#9003-07-0) Colour: White  
Odor: N/A N/A

### Section 4 : Normal Precautions

Ventilation	Not Necessary under normal conditions
Equipment	Not Necessary
Spills/Clean Up	Not Necessary

### Section 5: Fire and Explosion

Stability	Stable
Flash Point	N/A
Extinguishing Media	Foam, Carbon Dioxide, Water
Conditions to Avoid	None
Incompatible Materials	None

### Section 6: Other Precautions

Storage	No special precautions Necessary keep at normal room temperature
---------	---