

# C.H.A.R.L.I.E. Neonatal Resuscitation Simulator product profile



Product Name	<b>C.H.A.R.L.I.E. Neonatal Resuscitation Simulator</b>
Addresses Training Needs/Curricula	<ul style="list-style-type: none"><li>• Meets course curriculum for NRP® and Pediatric Advanced Life Support (PALS)</li><li>• Suitable for pre- and in-hospital training (EMS professionals and nurses)</li></ul>
Configuration	<ul style="list-style-type: none"><li>• Full-body manual simulator</li><li>• Multiple Packages for Optimal Training</li><li>• Heat-resistant (can be used in neonatal ICU unit and incubator)</li><li>• Features two birth anomalies for simulation</li><li>• SimVs offers an advanced simulation experience with C.H.A.R.L.I.E. NRP® Neonate with built-in scenario library and checklists</li></ul>
Dimension	<ul style="list-style-type: none"><li>• 19 1/2" (49.53 cm) from crown to heel</li><li>• Approximately 8 1/2 lbs (3.86 kg)</li></ul>
Product Code/SKU	LF01420 CHARLIE with Interactive ECG LF01421 CHARLIE

Medguard Professional Healthcare

Tel: 01 835 2378

Email: [orders@medguard.ie](mailto:orders@medguard.ie)

Fax: 01 969 5009

Web: [www.medguard.ie](http://www.medguard.ie)

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, Ireland



# Be READY

C.H.A.R.L.I.E. Neonatal Resuscitation Simulator


*Nasco*  
HEALTHCARE

**Ready to save lives.  
Right from the start.**

## **C.H.A.R.L.I.E. Neonatal Resuscitation Simulator**


Developed in conjunction with the latest Neonatal Resuscitation Program Guidelines (NRP®)

- C:** Compression and cardiac
- H:** Heat compatible: can be put in an incubator without melting
- A:** Airway, arterial access via the umbilicus
- R:** Resuscitation
- L:** Laryngeal mask capable for oxygen use or PPV
- I:** Intravenous, intraosseous, intubation
- E:** ECG



## **C.H.A.R.L.I.E. Neonatal Resuscitation Simulator:**

Full-term newborn simulator for training on the most common neonatal care scenarios: Enables nurses and EMS workers to practice the essentials of neonatal care.



This manikin is intended to be used as an instructional aide for teaching and demonstrating. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.