

Specification Sheet - FlowFlex Cov2 Rapid Antigen Test Nasal

Flowflex Cov2 Rapid Antigen Test Nasal



Relative Sensitivity	97.1% (83.8% -99.9%)
Relative Specificity	99.6% (97.7% -99.9%)
Accuracy	99.3% (97.5% -99.9%)
Test Format	Cassette
Specimen Format	Nasal Swab
Sample Volume	4 to 5 drops of the prepared sample
Test Time	15 Minutes
Shelf Life	24 Months Storage
Storage Temperature	2 – 30 ° C
Limit of Detection	1.6 * 10 ² TCID50 / ml



Overview:

- U.K. Government validated
- Full Home Test Approval in the UK and EU
- Detects New Virus Strains
- MHRA Registered
- Nasal only test
- 98.8% Accuracy
- CE approved and ISO certified
- 2 Year Shelf Life
- Results in 15 Minutes

Contents:

- 1 x Test Cassette
- 1 x Extraction Buffer Tube
- 1 x Disposable Swab
- 1 x Waste Bag
- 1 x Package Insert
- Instructions for use

How it works:

The FlowFlex Rapid Lateral Flow Test Kit uses the latest patented technology to detect proteins, ensuring that it can detect new strains of the COVID-19 virus. Using a simple nasal swab within 2cm of the nose makes it extremely quick & easy to administer an accurate test.

Use:

- Easy to use with fast results
- Affordable & accurate
- Suitable for health professionals & home testing
- Ideal for employees of all environments
- Suitable for regular, repeat testing

Order Number	Description
L031-118M5	FlowFlex Cov2 Rapid Antigen Test Nasal (1)

Please note: There are no returns on this item